

## EDTA Solution (0.2M in PBS)






Order No.	Size	Sterility
60-00030-10	10 ml	steril
60-00030-11	50 ml	steril
60-00030-12	100 ml	steril

### Beschreibung

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as an chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

### Additional Information

Interessengebiet	Zellbiologie
Sterilität	steril
ph-Bereich	pH 7.4
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.