

## ÜberStrainer Set 2

**Order No.:**43-70000-98



### Beschreibung

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 1 µm, 5 µm, 10 µm, 30 µm, 60 µm.

### Additional Information

Maschengröße	Multiple
Lieferumfang	10 Stück
Gewebematerial	PET (Polyethylenterephthalat)
Gehäusematerial	PP (Polypropylen)
Maschenart	gewoben
Sterilität	steril
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only. Not for use in diagnostic procedures.

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.