

Streck - Cell-Free DNA BCT® (10 ml, CE-IVD)



Order No.
Size

Order No.	Size
218996	6 tubes
218997	100 tubes
230244	1000 tubes

Beschreibung

A blood collection tube which stabilizes nucleated blood cells. The unique preservative prevents the release of genomic DNA, allowing isolation of high-quality cell-free DNA.

Convenient

- Reduces the need for immediate plasma preparation and CTC processing
- Allows for isolation of high-quality cell-free DNA for a wide range of downstream applications
- Allows for convenient sample collection, transport and storage

Stable

- cfDNA and gDNA are stable for up to 14 days at 6 °C to 37 °C
- CTCs are stable for up to 7 days at 15 °C to 30 °C
- Unfilled tubes have 24-month shelf life from the date of manufacture

Cell-Free DNA BCT is the trusted global leader in cell-stabilization

Additional Information

Vorgefüllt	Unfilled
Lieferzeit	1-2
Lagerbedingungen	Raumtemperatur

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

