

Syringe Strainer 40 µm (12 pcs., sterile)

Order No.:43-71040-51



Beschreibung

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

| | |
|-----------------------|--|
| Maschengröße | 40 µm |
| Variante | 12 pcs sterile |
| Lieferumfang | 12 pcs. |
| Durchflussregelung | Ja |
| Gehäusematerial | LD-PE (Polyethylen niedriger Dichte) |
| Maschenart | gewoben |
| Sterilität | steril |
| Lieferbedingungen | Raumtemperatur |
| Lagerbedingungen | Raumtemperatur |
| Regulierungserklärung | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.