

## Syringe Strainer 85 µm (12 pcs., sterile)

**Order No.:**43-71085-51



### Beschreibung

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Maschengröße	85 µm
Variante	12 pcs sterile
Lieferumfang	12 pcs.
Durchflussregelung	Ja
Gehäusematerial	LD-PE (Polyethylen niedriger Dichte)
Maschenart	gewoben
Sterilität	steril
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.