

Syringe Strainer 30 µm



Order No.	Sterility	Variant
43-71030-50	unsteril	12 pcs, non-sterile
43-71030-51	steril	12 pcs sterile

Beschreibung

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Maschengröße	30 µm
Lieferumfang	12 pcs.
Durchflussregelung	Ja
Gehäusematerial	LD-PE (Polyethylen niedriger Dichte)
Maschenart	gewoben
Sterilität	steril
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

