

Syringe Strainer 50 µm



| Order No. | Sterility | Variant |
|-------------|-----------|---------------------|
| 43-71050-50 | unsteril | 12 pcs, non-sterile |
| 43-71050-51 | steril | 12 pcs sterile |

Beschreibung

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

| | |
|-----------------------|--|
| Maschengröße | 50 µm |
| Lieferumfang | 12 pcs. |
| Durchflussregelung | Ja |
| Gehäusematerial | LD-PE (Polyethylen niedriger Dichte) |
| Maschenart | gewoben |
| Sterilität | steril |
| Lieferbedingungen | Raumtemperatur |
| Lagerbedingungen | Raumtemperatur |
| Regulierungserklärung | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

