

## Syringe Strainer 70 µm





Order No.	Sterility	Variant
43-71070-50	unsteril	12 pcs, non-sterile
43-71070-51	steril	12 pcs sterile

### Beschreibung

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

Maschengröße	70 µm
Lieferumfang	12 pcs.
Durchflussregelung	Ja
Gehäusematerial	LD-PE (Polyethylen niedriger Dichte)
Maschenart	gewoben
Sterilität	steril
Lieferbedingungen	Raumtemperatur
Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

