

Re-Strainer (In-Line-Strainer) 500 µm - Sterile

Order No.:43-75500-40



Beschreibung

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

Anwendung

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

Additional Information

Lieferzeit	1-2
Maschengröße	500 µm
Verpackung	4 devices per bag
Lieferumfang	20 pcs.

Farbe	neon pink
Durchflussregelung	ja
Mesh Fixation	Injected in housing for strong hold
Gewebematerial	PET (Polyethylenterephthalat)
Gehäusematerial	LD-PE (Polyethylen niedriger Dichte)
Sterilität	steril
Stabilität	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
Lieferbedingungen	Room Temperature
Lagerbedingungen	Raumtemperatur
Regulierungserklärung	For research use only
Rechtliche Hinweise	The product is for research and development only, not for diagnostic or therapeutic use

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.