

pluriSpin® Human T Cell Depletion



Order No.	Size
19-02002-31	2 ml
19-02002-32	10 ml

Beschreibung

pluriSpin® Human T Cell Depletion is designed to isolate pure, viable and untouched t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, while the isolated human t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Lieferzeit	1-2
Spezies	Mensch
Zelltypen	T-Zellen
Probenmaterial	Vollblut, PBMC, Buffy Coat, Nabelschnurblut, Knochenmark, primäre Zelllösung
Reinheit	>90%
Lagerbedingungen	4-8 °C
Regulierungserklärung	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

