

## Ly6G M-pluriBead® anti-ms



**Order No.:** 29-Ly6G0-20

### Description

Ly6G pluriBead® suspension for the positive cell isolation of Ly6G+ target cells from any murine sample material. (2ml)

### Additional Information

Size	2 ml
Species	Mouse (Murine)
Cell Types	Granulocytes
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.