

## CD19 M-pluriBead® anti-ms

**Order No.:**29-01900-20



### Description

CD19 pluriBead® suspension for the positive cell isolation of CD19+ target cells from any murine sample material. (2ml)

### Additional Information

|                      |                                                                               |
|----------------------|-------------------------------------------------------------------------------|
| Size                 | 2 ml                                                                          |
| Species              | Mouse (Murine)                                                                |
| Cell Types           | B cells                                                                       |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%                                                                          |
| Storage Condition    | 4-8 °C                                                                        |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.