

## pluriSpin® Human PBMC Depletion





Order No.	Size
19-08002-31	2 ml
19-08002-32	10 ml

### Description

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

### Additional Information

Species	Human
Cell Types	Non-PBMC cells
Storage Condition	Room Temperature
Regulatory Statement	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.