

ÜberStrainer Set 3

Order No.:43-70000-99



Description

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 100 µm, 200 µm, 300 µm, 400 µm, 500 µm.

Additional Information

| | |
|----------------------|--|
| Mesh Size | Multiple |
| Size | 10 pcs |
| Fabric material | PET (Polyethylenterephthalat) |
| Housing Material | PP (Polypropylene) |
| Mesh Type | woven |
| Sterility | Sterile |
| Shipping Condition | Room Temperature |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.