

## CD11b M-pluriBead® anti-ms

**Order No.:**29-011b0-20



### Description

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any murine sample material. (2ml)

### Additional Information

|                      |   |
|----------------------|---|
| Size                 | 2 ml  |
| Species              | Mouse (Murine)  |
| Cell Types           | Monocytes, Granulocytes, Macrophages, NK cells (subset)                       |
| Sample Material      | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Purity               | >90%  |
| Storage Condition    | 4-8 °C  |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.