

pluriSpin® Human CD45⁺ Cell Depletion



| Order No. | Size |
|-------------|-------|
| 19-09002-31 | 2 ml |
| 19-09002-32 | 10 ml |

Description

pluriSpin® Human CD45⁺ Cell Depletion is designed to isolate pure, viable and untouched CD45⁺-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45⁺ cells, while the isolated population free of CD45⁺ are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45⁺ cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

| | |
|----------------------|---|
| Species | Human |
| Cell Types | CD45 ⁻ Cell |
| Storage Condition | Room Temperature |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures., |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.