

## S-pluriBead® Mini Reagent Kit

**Order No.:**70-50010-11



### Descripción

Mini S-pluriBead Reagent Kit contains 100ml 10xWash buffer, stabilization buffer, incubation buffer, detachment activation buffer, 5x 1,8ml detachment buffer concentrate, 5 S-pluriStrainer, 5 Funnel, 5 connector rings

### Additional Information

Tamaño	5
Tiempo de entrega	1-2
Especies	Humano
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	Room Temperature; Buffer 4-8 °C; Buffer D -20°C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.