

CD4 M-pluriBead® anti-ms

Order No.:29-00400-20



Descripción

CD4 pluriBead® suspension for the positive cell isolation of CD4+ target cells from any murine sample material. (2ml)

Additional Information

Tamaño	2 ml
Especies	Ratón (Murine)
Tipos de células	Monocitos (subconjunto), Células T ayudantes
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.