

## CD8 M-pluriBead® anti-ms

**Order No.:**29-00800-20



### Descripción

CD8 pluriBead® suspension for the positive cell isolation of CD8+ target cells from any murine sample material. (2ml)

### Additional Information

Tamaño	2 ml
Especies	Ratón (Murine)
Tipos de células	Células NK (subconjunto), Células T citotóxicas
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.