

CD326 M-pluriBead® anti-ms

Order No.:29-32600-20



Descripción

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

Additional Information

Tamaño	2 ml
Especies	Ratón (Murine)
Tipos de células	Células tumorales en circulación, CTC (Células Tumorales en Circulación), Células epiteliales
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.