

## Ly6G M-pluriBead® anti-ms



**Order No.:**29-Ly6G0-20

### Descripción

Ly6G pluriBead® suspension for the positive cell isolation of Ly6G+ target cells from any murine sample material. (2ml)

### Additional Information

Tamaño	2 ml
Especies	Ratón (Murine)
Tipos de células	Granulocitos
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.