

## CD45 M-pluriBead® anti-hu

**Order No.:**19-04500-20



### Descripción

CD45 pluriBead® suspension for the positive cell isolation of CD45+ target cells from any human sample material. (2ml)

### Additional Information

|                             |   |
|-----------------------------|---|
| Tamaño                      | 2 ml  |
| Especies                    | Humano  |
| Tipos de células            | Leucocitos  |
| Material de muestra         | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Pureza                      | >90%  |
| Condición de almacenamiento | 4-8 °C  |
| Declaración reglamentaria   | For research use only. Not for use in diagnostic procedures,                  |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.