

## CD271 M-pluriBead® anti-hu

**Order No.:**19-27100-20



### Descripción

CD271 pluriBead® suspension for the positive cell isolation of CD271+ target cells from any human sample material. (2ml)

### Additional Information

Tamaño	2 ml
Especies	Humano
Tipos de células	Células mesenquimales, Células de Schwann, Células del tallo/estromal
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.