

CD11b M-pluriBead® anti-hu

Order No.:19-011b0-20



Descripción

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any human sample material. (2ml)

Additional Information

Tamaño	2 ml
Especies	Humano
Tipos de células	Monocitos, Granulocitos, Macrófagos, Células NK (subconjunto)
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.