

CD11b S-pluriBead® anti-ms

Order No.:29-011b0-10



Descripción

CD11b pluriBead® suspension for the positive cell isolation of CD11b+ target cells from any murine sample material. (2ml)

Additional Information

Tamaño	2 ml
Especies	Ratón (Murine)
Tipos de células	Monocitos, Granulocitos, Macrófagos, Células NK (subconjunto)
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.