

CD55 S-pluriBead® anti-hu

Order No.:19-05500-10



Descripción

CD55 pluriBead® suspension for the positive cell isolation of CD55+ target cells from any human sample material. (2ml)

Additional Information

| | |
|-----------------------------|---|
| Tamaño | 2 ml |
| Especies | Humano |
| Tipos de células | Células hematopoiéticas, Células no hematopoyéticas |
| Material de muestra | Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution |
| Pureza | >90% |
| Condición de almacenamiento | 4-8 °C |
| Declaración reglamentaria | For research use only. Not for use in diagnostic procedures, |

Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.