

## 10x PBS Solution





Order No.	Size	Sterility
60-00010-10	100 ml	Estéril
60-00010-11	500 ml	Estéril

### Descripción

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without**:

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (1368,9 mMol/L)
- Na<sub>2</sub> HPO<sub>4</sub> • 2H<sub>2</sub> O (100mMol/L)
- KH<sub>2</sub> PO<sub>4</sub> (19,8 mMol/L)
- KCl (26,8 mMol/L)

### Additional Information

Tiempo de entrega	1-2
Área de interés	Biología celular
Esterilidad	Estéril
Rango de pH	pH 6.8
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación

Declaración reglamentaria

For research use only. Not for use in diagnostic procedures,

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.