

## 10% BSA Stock Solution in 1x PBS






Order No.	Size	Sterility
60-00020-10	50 ml	Estéril
60-00020-11	250 ml	Estéril
60-00020-12	500 ml	Estéril

### Descripción

10% BSA stock solution in 1x PBS can used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

### Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

### Additional Information

Tiempo de entrega	1-2
Área de interés	Biología celular
Esterilidad	Estéril
Rango de pH	pH 7.4
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

