


## EDTA Solution (0.2M in PBS)






Order No.	Size	Sterility
60-00030-10	10 ml	Estéril
60-00030-11	50 ml	Estéril
60-00030-12	100 ml	Estéril

### Descripción

EDTA Stock Solution is a sterile, pH adjusted 100x concentrate. Dilute the concentrate to achieve the working reagent. EDTA binds calcium and other metals and is usable as a chelating agent for preserving blood sample. Moreover EDTA can be added to buffers for cell washing steps or in running buffers for flow cytometry or automated cell separation methods.

### Additional Information

Área de interés	Biología celular
Esterilidad	Estéril
Rango de pH	pH 7.4
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.