

pluriSpin® Human CD45⁺ Cell Depletion (10 ml)

Order No.:19-09002-32



Descripción

pluriSpin® Human CD45⁺ Cell Depletion is designed to isolate pure, viable and untouched CD45⁺-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels CD45⁺ cells, while the isolated population free of CD45⁺ are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD45⁺ cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Tamaño	10 ml
Especies	Humano
Tipos de células	CD45- Célula
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures.,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.