

## pluriSpin® Human CD4<sup>+</sup> T Cell Enrichment





Order No.	Size
19-03001-31	2 ml
19-03001-32	10 ml

### Descripción

pluriSpin® Human CD4<sup>+</sup> T Cell Enrichment is designed to isolate pure, viable and untouched CD4<sup>+</sup> T cells from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels Monocytes, cytotoxic T cells, B cells, NK cells and granulocytes, while the negatively isolated human CD4<sup>+</sup> T cells are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified CD4<sup>+</sup> T cells are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Especies	Humano
Tipos de células	Células T ayudantes
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.