

## pluriSpin® Human PBMC Depletion (10 ml)

Order No.:19-08002-32



### Descripción

pluriSpin® Human PBMC Depletion is designed to the Peripheral Blood Mononuclea Cell (PBMC T cells, B cells, NK cells, Monocytes) population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels PBMC. When centrifuged over a density medium the unwanted cells pellet along with the RBCs.

### Additional Information

Tamaño	10 ml
Especies	Humano
Tipos de células	Células no PBMC
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.