

## Centrifuge Tube 50 ml (sterile, 300 pcs.)



**Order No.:**05-00001-01

### Descripción

50 ml tubes are suitable for a wide range of laboratory applications such as sample preparation or sample storage. The tubes are ideal for the use with our **pluriStrainers** and guarantee a perfect fit with the Connector Ring. Moreover the tubes can be used as container for **density gradient centrifugation** or the enrichment of primary cells with **pluriBead**.

Screw tube, 50 ml, 114 x 28 mm, pointed bottom, PP, with printed writing field and scaling, with mounted red cap, sterile and pyrogen-free / endotoxin-free, 25 pcs./bag

### Aplicación

- Mixing
- Incubation
- Centrifugation
- Storage
- Collection vessel for filtration and purification of liquids

### Additional Information

Tamaño	300 pcs.
Pre-llenado	Unfilled
Tiempo de entrega	1-2
Variación	20 pcs, sterile
Embalaje	25 pcs. / Bag
Condición de almacenamiento	La temperatura de la habitación
Condición de envío	La temperatura de la habitación

Declaración reglamentaria

For research use only. Not for use in diagnostic procedures.

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.