

Bundle pluriStrainer® Mini 5 µm and Flow Cytometry Tubes



	Order No.	Size	Sterility	Packaging
<input type="checkbox"/>	43-10005-46	25 pcs.	Estéril	25 pcs. (25 pcs. / Bag)
<input type="checkbox"/>	43-10005-56	50 pcs.	Estéril	50 pcs. (25 pcs. / Bag)
<input type="checkbox"/>	43-10005-66	100 pcs.	Estéril	100 pcs. (25 pcs. / Bag)
<input type="checkbox"/>	43-10005-76	500 pcs.	Estéril	500 pcs. (25 pcs. / Bag)

Descripción

Bundle includes sterile pluriStrainers Mini with 5 µm mesh and sterile 5 ml polystyrene (PS) tubes (12x75mm). The combination of a flow cytometry tube and cell strainer offers a convenient way to prepare laboratory samples for flow cytometry analysis. A 5 µm polyethyleneterephthalat (PET) mesh is incorporated into the housing of the pluriStrainer Mini 5 µm. The polystyrene round-bottom tubes are closed with a cap. Tubes and the mini cell strainers will be delivered in separate bags with 25 pcs. each.

Aplicación

- Obtaining a real single cell suspension after digestion of mammary tissue and organoids
- Preparation of a single cell suspension of blood cells from bone marrow, pancreas, thymus, lymph nodes and others
- Dissociation of cells from other primary tissue
- Preparation of real single cell suspension for flow cytometry (FACS™)

Additional Information

El tamaño de la malla	5 µm
Embalaje	50 pcs. (25 pcs. / Bag)
Color	amarillo
Control de flujo	No
Mesh Fixation	Injected in housing for strong hold
Material de la tela	PET (Polietilentereftalato)
Material de la carcasa	PP (Polipropileno)

Encaja en	1.5 ml reaction tube, 2.0 ml reaction tube, FACS™ tube, Cryo vials 5 ml / 10 ml, reaction tube with lid, 5 ml / 10 ml reaction tube with screw cap, 14- 15 ml conical centrifuge tube, 24 well plate, 48 well plate
Manejo	Gripping surface enables easy and aseptic handling.
Esterilidad	Estéril
Comparable con	Cell Strainer Snap Cap™, Corning Cell Strainer, EASYstrainer™, Cell Sieve, Falcon™ Cell Strainer, Reversible™ Strainer, MACS® SmartStrainer, MACS® Pre-Separation Filters
Estabilidad	RLT® Buffer, TRIzol®, isopropanol, organic solvents
Centrifugable	Sí
Disociación de tejidos	Sí
Desinfectable	Sí
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Información legal	Buffer RLT® is a trademark of Qiagen, TRIzol® is a trademark of Molecular Research Center Inc., FACS™ is a trademark of BD Biosciences, Corning is a registered trademark of Corning Inc., EASYstrainer™ is a registered trademark of Greiner Bio One International GmbH, Cell Strainer Snap Cap™ and Falcon™ Cell Strainers is a registered trademark of Corning Inc., MACS® SmartStrainer is a registered trademark of Miltenyi Biotec GmbH

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.