

## ÜberStrainer Set 1

**Order No.:**43-70000-97



### Descripción

The ÜberStrainer is more than a cell strainer - ÜberStrainer is a sterile Sample Preparation Device to separate, isolate, concentrate and process particle out of complex and/ or diluted suspensions. The air tight screw cap with a Luer Lock connector makes it possible to filtrate an unlimited amount of liquid. The strainer can be easily separated from the rest of the device, fits into a 50 ml tube and into the 20 mm tissue plate cavity.

This set includes 2 pieces of 20 µm, 60 µm, 100 µm, 300 µm, 500 µm.

### Additional Information

El tamaño de la malla	Multiple
Tamaño	10 pcs
Material de la tela	PET (Polietilentereftalato)
Material de la carcasa	PP (Polipropileno)
Mesh Type	woven
Esterilidad	Estéril
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures.

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.