

## pluriMag Rack I (1.5 / 2 ml tubes)

**Order No.:**50-01010-01



### Descripción

pluriMag Rack is a magnetic holder for small-scale protein and nucleic acids purification and cell enrichment or depletion with magnetic beads. pluriMag Rack can be used for sample volumes of up to 5 ml. The holder make an excellent complement to the following pluriMag products (containing magnetic beads) from pluriSelect. pluriMag Rack works also with all standard magnetic micro beads.

It is optimized for efficient magnetic separation of all types of magnetic (micro) beads (range 1–8  $\mu\text{m}$  in diameter) from the supernatant (10–2,000  $\mu\text{l}$ ) in seconds, allowing for easy pipette removal without sample loss. The pluriMag Rack combines a strong magnetic field for efficient separations in standard microcentrifuge tubes.

- Best working volume: 10–2,000  $\mu\text{l}$  (1ml -50ml)
- Holds up to 10 standard 1.5–2 ml microcentrifuge tubes in or 4 flow cytometry tubes
- Good control and visibility of your sample

### Additional Information

Tamaño	1 device
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures.,

### Warning and Limitations

This product is for research and development only, not for diagnostic or theurapeutic use.