

Syringe Strainer 5 µm (12 pcs., non-sterile)

Order No.:43-71005-50



Descripción

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

Tiempo de entrega	1-2
El tamaño de la malla	5 µm
Variación	12 pcs, non-sterile
Tamaño	12 pcs.
Control de flujo	Sí
Material de la carcasa	LD-PE (Polietileno de baja densidad)
Mesh Type	woven
Esterilidad	No estéril
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures.

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.