

Syringe Strainer 100 µm (12 pcs., non-sterile)

Order No.:43-71100-50



Descripción

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

Additional Information

| | |
|-----------------------------|--|
| El tamaño de la malla | 100 µm |
| Variación | 12 pcs, non-sterile |
| Tamaño | 12 pcs. |
| Control de flujo | Sí |
| Material de la carcasa | LD-PE (Polietileno de baja densidad) |
| Mesh Type | woven |
| Esterilidad | No estéril |
| Condición de envío | La temperatura de la habitación |
| Condición de almacenamiento | La temperatura de la habitación |
| Declaración reglamentaria | For research use only. Not for use in diagnostic procedures. |

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.