

## Syringe Strainer 5 µm (12 pcs., sterile)



**Order No.:**43-71005-51

### Descripción

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

|                             |  |
|-----------------------------|--|
| Tiempo de entrega           | 1-2  |
| El tamaño de la malla       | 5 µm   |
| Variación                   | 12 pcs sterile   |
| Tamaño                      | 12 pcs.  |
| Control de flujo            | Sí   |
| Material de la carcasa      | LD-PE (Polietileno de baja densidad)                         |
| Mesh Type                   | woven  |
| Esterilidad                 | Estéril  |
| Condición de envío          | La temperatura de la habitación                              |
| Condición de almacenamiento | La temperatura de la habitación                              |
| Declaración reglamentaria   | For research use only. Not for use in diagnostic procedures. |

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.