

## Syringe Strainer 10 µm





Order No.	Sterility	Variant
43-71010-50	No estéril	12 pcs, non-sterile
43-71010-51	Estéril	12 pcs sterile

### Descripción

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

El tamaño de la malla	10 µm
Tamaño	12 pcs.
Control de flujo	Sí
Material de la carcasa	LD-PE (Polietileno de baja densidad)
Mesh Type	woven
Esterilidad	Estéril
Condición de envío	La temperatura de la habitación
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures.

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

