

## 1x PBS Solution (100 ml)

**Order No.:**60-00000-01



### Descripción

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub> HPO<sub>4</sub> • 2H<sub>2</sub> O (10mMol/L)
- KH<sub>2</sub> PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Additional Information

Tamaño	100 ml
Tiempo de entrega	1-2
Área de interés	Biología celular
Esterilidad	Estéril
Rango de pH	pH 7.4
Condición de envío	La temperatura de la habitación

Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.