

## Re-Strainer (In-Line-Strainer) 1 µm - Sterile

**Order No.:**43-75001-40



### Descripción

The Re-Strainer is a filtration device for high volumes of liquid to exclude or concentrate particles from liquids with the possibility of particle recovery. The 6 ml recovery reservoir can be accessed by unscrewing the cap. The Re-Strainer has two female luer-lock-ports.

These ports can be used to connect several Re-Strainer to a filter cascade (male to male luer-lock connector) or integrate the Re-Strainer as an in-line-strainer via luer-lock tube adaptor.

The mesh is an integrated part of the housing to give the device maximum mechanical stability.

### Aplicación

- Filtration of large sample volume
- Cascade filtration
- Size-fractionation of particles
- Removal of smaller impurities
- Concentration of rare particles from large volumes
- Cartridge for affinity chromatography gels

### Additional Information

Tiempo de entrega	1-2
El tamaño de la malla	1 µm
Variación	20 pcs, sterile
Embalaje	4 devices per bag

Tamaño	20 pcs.
Color	rojo
Control de flujo	sí
Mesh Fixation	Injected in housing for strong hold
Material de la tela	PET (Polietilentereftalato)
Material de la carcasa	LD-PE (Polietileno de baja densidad)
Esterilidad	Estéril
Estabilidad	Good resistance to non-oxidizing acids and bases, fats, most organic solvents
Condición de envío	Room Temperature
Condición de almacenamiento	La temperatura de la habitación
Declaración reglamentaria	For research use only
Información legal	<b>The product is for research and development only, not for diagnostic or therapeutic use</b>

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.