

pluriSpin® Human T Cell Depletion



Order No.	Size
19-02002-31	2 ml
19-02002-32	10 ml

Descripción

pluriSpin® Human T Cell Depletion is designed to isolate pure, viable and untouched t cells-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels T cells, while the isolated human t cells-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified t cells-free population are present as a highly enriched population at the interface between the plasma and the density medium.

Additional Information

Tiempo de entrega	1-2
Especies	Humano
Tipos de células	Células T
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

