

## pluriSpin® Human Monocyte Depletion (2 ml)

**Order No.:**19-01002-31



### Descripción

pluriSpin® Human Monocyte Depletion is designed to isolate pure, viable and untouched monocytes-free population from whole blood, buffy coat, cord blood or similar sample material by negative isolation. The product labels monocytes, while the isolated human monocytes-free population are left in the sample. When centrifuged over a density medium the unwanted cells pellet along with the RBCs. The purified monocytes-free population are present as a highly enriched population at the interface between the plasma and the density medium.

### Additional Information

Tamaño	2 ml
Especies	Humano
Tipos de células	Monocitos
Material de muestra	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureza	>90%
Condición de almacenamiento	4-8 °C
Declaración reglamentaria	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.

