


## 10% BSA stock solution in 1x PBS (500 ml)

**Order No.:**60-00020-12



### Description

10% BSA stock solution in 1x PBS can be used for cell culture, cell separation wash buffer, blocking steps in Western blot, ELISA, IHC and nucleic acid detection methods.

### Features of Blocker BSA

- **Ready to use** - No need to dissolve powder
- **Pure** - High-quality bovine serum albumin as 10% solution provides fewer chances of cross-reaction with assay components than other components
- **Convenient** - Optimal blocking results for specific applications

### Additional Information

|                      |  |
|----------------------|--|
| Size                 | 500 ml   |
| Delivery Time (days) | 1-2  |
| Area of interest     | Cell Biology   |
| Sterility            | Sterile  |
| pH Range             | pH 7.4   |
| Shipping Condition   | Room Temperature   |
| Storage Condition    | 4-8 °C   |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures, |

## **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.