

## CD55 M-pluriBead® anti-hu

**Order No.:**19-05500-20



### Description

CD55 pluriBead® suspension for the positive cell isolation of CD55+ target cells from any human sample material. (2ml)

### Additional Information

Size	2 ml
Species	Human
Cell Types	Hematopoietic cells, Non-hematopoietic cells
Sample Material	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Purity	>90%
Storage Condition	4-8 °C
Regulatory Statement	For research use only. Not for use in diagnostic procedures,

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.