

## PLT Spin Medium (100 ml)

**Order No.:** 60-00094-10



### Description

PLT Spin Medium is a ready-to-use, sterile density medium for the isolation of platelets in high yield from fresh human peripheral blood through single step centrifugation.

### Additional Information

Size	100 ml
Delivery Time (days)	1-2
Storage Condition	Room temperature, protect from light
Regulatory Statement	For research use only. Not for use in diagnostic procedures.
pH Range	6.8 – 7.4

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.