

## Syringe Strainer 20 µm (12 pcs., sterile)

**Order No.:**43-71020-51



### Description

Sieving or filtration device for large quantities to exclude or concentrate particles from liquids. The strainer can be connected to a syringe or a tube via the screw cap with a Luer-Lock adaptor which works under positive or negative pressure.

### Additional Information

|                      |  |
|----------------------|--|
| Mesh Size            | 20 µm  |
| Variation            | 12 pcs sterile   |
| Size                 | 12 pcs.  |
| Flow Control         | Yes  |
| Housing Material     | LD-PE (Low Density Polyethylen)                              |
| Mesh Type            | woven  |
| Sterility            | Sterile  |
| Shipping Condition   | Room Temperature   |
| Storage Condition    | Room Temperature   |
| Regulatory Statement | For research use only. Not for use in diagnostic procedures. |

### Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.