

## 1x PBS Solution (500 ml)

**Order No.:**60-00000-10



### Description

PBS (phosphate buffered saline) is a balanced salt solution used for a variety of cell culture applications, such as washing cells before dissociation, transporting cells or tissue, diluting cells for counting and preparing reagents. PBS is formulated without calcium and magnesium for rinsing chelators from the culture before cell dissociation.

This PBS is manufactured **without:**

- Calcium
- Magnesium
- Phenol Red
- Protease free

### Components

- NaCl (136,89 mMol/L)
- Na<sub>2</sub> HPO<sub>4</sub> • 2H<sub>2</sub> O (10mMol/L)
- KH<sub>2</sub> PO<sub>4</sub> (1,98 mMol/L)
- KCl (2,68 mMol/L)

### Informations complémentaires

Taille	500 ml
Délai de livraison	1-2
Domaine d'intérêt	Biologie cellulaire
Stérilité	Stérile
Gamme de ph	pH 7.4
Conditions d'expédition	Température ambiante

Les conditions de stockage	Température ambiante
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

### **Warning and Limitations**

This product is for research and development only, not for diagnostic or therapeutic use.