

CD326 M-pluriBead® anti-ms

Order No.:29-32600-20



Description

CD326 pluriBead® suspension for the positive cell isolation of CD326+ target cells from any murine sample material. (2ml)

Informations complémentaires

Taille	2 ml
Espèces	Souris (Murine)
Types de cellules	Circulation des cellules tumorales, CTC (cellules tumorales circulantes), Cellules épithéliales
Échantillon de matériel	Whole Blood, PBMC, Buffy Coat, Cord Blood, Bone Marrow, Primary Cell Solution
Pureté	>90%
Les conditions de stockage	4-8 °C
Déclaration réglementaire	For research use only. Not for use in diagnostic procedures,

Warning and Limitations

This product is for research and development only, not for diagnostic or therapeutic use.